


# COVID-19 Ag Detection Kit (Immunofluorescence-Based)

## ● WHEN TO USE ANTIGEN TEST?




## ● INTENDED USE

COVID-19 Ag Detection Kit (Immunofluorescence-Based) is an in vitro rapid test for the detection of SARS-CoV-2 antigen in direct nasal swabs from individuals suspected of COVID-19 within the first seven days of symptom onset.

COVID-19 Ag Detection Kit (Immunofluorescence-Based) is intended for use by medical professionals or trained operators who are proficient in performing rapid lateral flow tests.

## ● ASSAY PROCEDURES:




## ● INTERPRETATION OF TEST RESULT

### Dry Immunofluorescence Analyzer

Positive			Negative
Weak	Medium	Strong	
$0.25 \geq T/C > 0.15$	$2.5 \geq T/C > 0.25$	$T/C > 2.5$	$\leq 0.15$
+	++	+++	-

### Blue LED Light



## ● PERFORMANCE:

### Clinical Performance

Reagents		RT-PCR		Total
		Positive	Negative	
COVID-19 Ag Detection Kit	Positive	74	1	75
	Negative	5	131	136
Total		79	132	211

Relative Sensitivity: 93.67%, 95% CI: 88.41%~98.93%;  
Relative Specificity: 99.24%, 95%CI: 95.38%~99.31%.

### Analytical Sensitivity

Concentration TCID <sub>50</sub> /ml	Concentration copies/ml	Number Positive/Total	% Detected
93.3	$5.67 \times 10^5$	20/20	100

The analytical sensitivity of COVID-19 Ag Detection Kit (Immunofluorescence-Based) is 93.3 TCID<sub>50</sub>/ml.

Further information:  
[www.torontobioscience.com](http://www.torontobioscience.com) and  
[sales@torontobioscience.com](mailto:sales@torontobioscience.com)

# COVID-19 Ag Detection Kit (Immunofluorescence-Based)



Further information:

[www.torontobioscience.com](http://www.torontobioscience.com) and [sales@torontobioscience.com](mailto:sales@torontobioscience.com)


**IMPORTANT:** Read the instruction manual before use. **For the best performance, direct nasal swabs should be tested preferably as soon as possible after collection.**

## PROCEDURE CARD

### ASSAY PROCEDURE


#### Specimen Preparation Procedure: (Nasal swab)

- Transfer 20 drops (~ 500  $\mu\text{L}$ ) of sample lysis buffer into an extraction tube using a transfer pipet.
- Insert the swab into an extraction tube. While squeezing the buffer tube, stir the swab more than 5 times and wait for 1 minute.
- Squeeze the wall of the tube to extract the liquid from the swab.
- Press the nozzle cap tightly onto the tube.



#### Analysis of Specimen:

- Apply 2 drops (~80  $\mu\text{L}$ ) of the extracted specimen to the specimen well of the test strip.
- Read the test result in 20 minutes.



### RESULT INTERPRETATION

#### Dry Immunofluorescence Analyzer Run Test with FIC-H1W

- Turn on the machine
- Select "QUICK TEST"
- Input patient information, including Name, Age, Sex.
- Insert the prepared test strip into the machine.
- Click "QUICK TEST"
- The test result (RS) will be displayed on the screen within 5 seconds.


Positive			Negative
Weak	Medium	Strong	
$0.25 \geq T/C > 0.15$	$2.5 \geq T/C > 0.25$	$T/C > 2.5$	$\leq 0.15$
+	++	+++	-

**Note:** If the result is invalid, the test should be repeated.

#### Blue LED Light

If Portable fluorescent immunoassay analyzer is not available, a blue LED light also can be used to read the test result.

Wears glasses when reading results using blue LED light.



Positive: the presence of two lines as the control line (C) and the test line (T) in the result window.

Negative: the presence of a single line as the control line (C) in the result window.

Invalid: if the control line (C) is not visible within the result window after performing the test, the result is invalid.